



High Resolution
Fibre Endoscopy



ENT Scopes

A diverse range of Conventional
and Portable fibreoptic scopes
for high quality ENT care

PENTAX

High Resolution Fibre ENT Scopes



A diverse and unique choice of Conventional and Portable ENT fibrescopes for a wide range of applications. Superior quality with enhanced observation and handling capabilities delivering reliable and safe results every time.

Outstanding Performance Characteristics



Maximum Working Comfort

- 130° tip deflection, a tight bending radius and up to 100° angle of view for easy manoeuvrability and more precise and accurate examinations.
- Ergonomically designed, lightweight control body enhancing user performance.
- Specially rounded distal tip eases introduction and minimises patient discomfort.
- PENTAX LSV Solutions (Light Source Versatility). Three Light Source options; The Cordless battery powered LSV module or self contained A/C adapter Light Source module (FNL-10RBS only). Finally all models can be connected directly to a conventional light source.

Unique Hygiene

- Repositioned channel* inlet away from eyepiece improves user safety and reduces risk of contamination.
- Fully immersible with new PENTAX Safety Seal for improved reprocessing performance.

* Only in FNL-15RP3

Superb Image Quality

- PENTAX Superfine Fibre Technology delivering bright, clear and high resolution images for superior observation capabilities.
- Upgradeable to video image through connection of ocular mounted camera system.

PENTAX Quality & Service Programme

- The PENTAX brand continues to be recognised for delivering cost effective products of exceptional quality and durability.
- Industry leading service plans, consultancy, financing, technical service and training.

Diagnostic/Therapeutic Fibre ENT Scopes

| Type | Ultra-Slim | Routine | Portable | Therapeutic |
|----------------------------------|------------|-----------|-----------|-------------|
| Model | FNL-7RP3 | FNL-10RP3 | FNL-10RBS | FNL-15RP3 |
| Angle of view (°) | 75 | | 85 | 75 |
| Focal range (mm) | 3-50 | | | |
| Tip deflection (°) up/down | 130/130 | | | |
| Insertion tube diameter (mm) | 2.4 | | 3.4 | 4.8 |
| Distal end diameter (mm) | 2.4 | | 3.5 | 4.9 |
| Instrument channel diameter (mm) | - | - | - | 2.2 |
| Working length (mm) | 300 | | | |

PENTAX Europe GmbH
LIFE CARE
Julius-Vosseler-Straße 104
22527 Hamburg
Deutschland
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: medical@pentax.de

PENTAX U.K. Limited
LIFE CARE
Pentax House
Heron Drive, Langley
Slough SL3 8PN
United Kingdom
Tel.: +44 17 53 / 79 27 23
Fax: +44 17 53 / 79 27 94
E-mail: lifecare.sales@pentax.co.uk

PENTAX France S.A.S.
LIFE CARE
112 quai de Bezens
P. B. 204
95106 Argenteuil
France
Tel.: +33 1 / 30 25 75 75
Fax: +33 1 / 30 25 75 76
E-mail: medical@pentax.fr

PENTAX Nederland B.V.
LIFE CARE
Lage Mosten 35
4822 NK Breda
Netherlands
Tel.: +31 76 / 5 31 30 31
Fax: +31 76 / 5 31 30 00
E-mail: lifecare@pentax.nl

PENTAX Italia S.r.l.
LIFE CARE
Via Dione Cassio, 15
20138 Milano
Italy
Tel.: +39 / 02 50 99 58 1
Fax: +39 / 02 50 99 58 60
E-mail: marketing.lifecare@pentaxitalia.it

S.I.M. - Sistemas
Integrales de Medicina, S.A.
LIFE CARE
Sistema Solar 25
28830 San Fernando de Henares
Madrid · Spain
Tel.: +34 91 / 301 62 40
Fax: +34 91 / 751 31 15
E-mail: sim@simmedica.com

HOYA Corporation
PENTAX Tokyo Office
2-36-9, Maeno-cho
Itabashi-Ku
174-8639 Tokyo
Japan
Tel.: +81 33 / 9 60 51 55
Fax: +81 35 / 3 92 67 24

www.pentaxeurope.com



0197 HOYA Corporation · ISO 9001/EN ISO 13485 certified

Your PENTAX sales representative:

PENTAX