

Accessories Programme



Periphery

Everything you need to customise your PENTAX endoscope system.

PENTAX



The PENTAX accessories programme of peripheral equipment for your endoscope

Greater degree of compatibility

Finely tuned system components are essential for safe and successful endoscopy. That is why we aim to provide the highest possible quality in every detail of the peripheral equipment for PENTAX endoscopes and endoscopy systems. The following pages describe a carefully compiled accessory programme for customising and installing additional equipment for your PENTAX endoscopy system.



How to choose the correct accessories:

Before you decide on new accessories, we recommend that you read the following three-point checklist. Our PENTAX specialists will be pleased to advise you if you have any questions regarding your endoscope and suitable accessories for it.

- Please check that your endoscope is designed for/compatible with the desired peripheral equipment. Important relevant information can be found on your endoscope's or processor's identification plate.
- When choosing new instruments, please be sure to note the diameter of the channel (see identification plate) and the available working length of your endoscope.
- All PENTAX peripheral equipment in this catalogue has been tested in conjunction with PENTAX endoscopes for functionality and compatibility.

Suction pumps and accessories	4
Video Classic carts and accessories	5
Swivelling arm	6
Documentation equipment and accessories	6-7
Cable	
● Video Cable.....	8
● Remote control cable.....	8
● Various Cable	8
● Multi-Video-Hybrid-Cable	9
● Cable for High Frequency Instruments	9
Adapter	
● Olympus	10-13
● Storz.....	14
● Wolf	14
● Rigid Endoscopes	15
● Eyepiece Accessory	15
Adapter PENTAX	
● PENTAX Adapter for Processor, Videomodule, Lightsources.....	16
Accessories	
● for Bedsidescopes	17
● for Halogen light sources.....	18
● for Lightsources/Processors	18
● for Luer Lock fitting	19
● for Gasinsuflation.....	19

Suction pumps and accessories

Senator-Pump



Code	Article	Description
300065	Senator 30	Senator 30 suction pump incl. 2 litre container, troughput: 30 litre/min
300062	2L Container	2L container, autoclavable for Senator 30 pump
300063	3L Container	3L container, autoclavable for Senator 30 pump
300064	5L Container	5L container, autoclavable for Senator 30 pump
300135	Cover	Cover with handle for 2L,3L and 5L container/Senator pump
300136	Filter	Filter for Senator 30 pump
300071	Holder	Holder for Senator 30 Pump

Medap-Pump



Code	Article	Description
300066	Medap	MEDAP suction pump type P7050
300086	1L Container	1L container, autoclavable for MEDAP pump




Video Classic carts and accessories

Video Classic carts and accessories



Code	Article	Description
400700	ITDMITRAFO	ITD Classic-cart PENTAX with transformator
400701	ITDOHNETRAFO	ITD Classic-cart PENTAX without transformator
400702	ITDCONFOCAL	ITD Classic-cart for PENTAX Endomicroscopy-system
400702	ITDTRAFOISO	ITD Classic-cart PENTAX with transformator and isoguard

Technical data	Classic-cart
Weight (kg)	70
Dimensions (H x L x W) mm	1,277 x 614 x 625
Total load (kg)	150





Swivelling arm

Swivelling arm for ITD trolley

Code	Article	Description
400517	TRHA-01	Swivelling arm for ITD trolley for 15" and 19" Monitor

Sony Base

Code	Article	Description
403022	SU560	Base for Sony LMD-1950MD and LMD-2140MD

Wall holder

Code	Article	Description
400523	WHNO02	Wall holder

Documentation equipment

Color Videoprinter

Code	Article	Description
400125	UP-21MD	Color videoprinter, Sony

DVD-Recorder

Code	Article	Description
403396	LQ-MD800E	Panasonic DVD-RAM /R Recorder

MediCap-System

Code	Article	Description
400008	USB100	Picture-Storage-unit incl. 1 USB-Stick (256MB)
400009	USB200	Picture-and Video-Storage-unit incl. 1 USB-Stick (512MB)

HD-XDCAM System (for storing video sequences in HD quality)

Code	Article	Description
-	-	SONY PDW70MD (on request)

Documentation equipment and accessories

Documentation equipment and accessories

Code	Article	Description
403361	CK-700	Paper for Mitsubishi CP-700/710E (PK-700S and PK-700L)
403362	PK-700L	Ink-Sheet For CP-700E/710E (120 prints large format)
403363	PK-700S	Ink-Sheet For CP-700E/710E (200 prints normal format)
403371	CK-900L	Printing paper and Ink-Sheet For CP-900 (large)
403370	CK-900S	Printing paper and Ink-Sheet For CP-900 (short)
400098	UPC-2010	Printing paper and Ink-Sheet for UP-2000-series, 200 sheet
400105	UPC-1010E	Printing paper and Ink-Sheet for UP-1850/1800/1200, 100 sheet
400102	UPC-3010P	Printing paper and Ink-Sheet for UP-3000-series, 80 sheet (large format)
400128	UPC-21L	Printing paper and Ink-Sheet for UP-21MD, 200 sheet (short format)
400127	UPC-21S	Color-print-package for UP-21MD, 240 prints possible
403355	CK-50S	Printing Paper for Mitsubishi CP-50E, 50 sheet
403185	Switch	Switch Y/C 4 input to 1 output
400107	FS-20	Sony foot-operated switch to activate video printers

Cable

Videocable

Code	Article	Description	Length
400560	DVI-D to DVI-D Cable	Connection between EPK-i and monitor	1.8 m
400561	DVI-D to DVI-D Cable	Connection between EPK-i and monitor	3 m
400565	DVI Fiber Optic Cable	Connection between EPK-i and monitor	10 m
400562	DVI Fiber Optic Cable	Connection between EPK-i and monitor	20 m
81532	OS-A25	D-SUB RGBS Cable for EPK-700/1000/i, EPM-1000/3500	1.3 m
403180	RGB 9 pol. 5 m	RGB-Cable 9 polig D-SUB/BNC for EPK-700/1000/i, EPM-1000	5 m
403181	RGB 9 pol. 10 m	RGB-Cable 9 polig D-SUB/BNC for EPK-700/1000/i, EPM-1000	10 m
403158	RGB 1,2 m	RGB-Cable BNC/Scart	1.2 m
403156	RGB 10 m BNC Scart	RGB-Cable BNC/Scart	10 m
403150	RGB 5 m	RGB-Cable BNC/BNC	5 m
403154	RGB 10 m	RGB-Cable BNC/BNC	10 m
403176	S-VHS 2 m	S-VHS Cable	2 m
403175	S-VHS 5 m	S-VHS Cable	5 m
403172	S-VHS 10 m	S-VHS Cable	10 m

Remote Control Cables

Code	Article	Description	Length
86364	OS-A58	Remote control cable for processors EPK-1000/100p/i	-
85998	OS-A48	Remote control cable for processor	-
85997	OS-A47	Control cable for printer (out)	-
403147	Remote Control Cable	Remote control cable for video printer processor 3.5 mm plug; double-sided; stereo for EPK-100 & EPM-3500	-
403146	Remote Control Cable	Remote control cable for video recorder Panasonic AG-MD800, 2.5 mm mini-plug	1.5 m

Various Cable

Code	Article	Description	Length
86389	OL-Z3	Earthing box EPK-1000	1,3 m
85995	OS-A45	Extension cable for keyboard	-
400800	POAG60	POAG cable	0,6 m
400801	POAG300	POAG cable	3 m
400802	POAG500	POAG cable	5 m



Multi-Video Hybrid Cable for EPK-700/1000



Code	Article	Description
403144	Multi-Video Hybrid Cable 10 m for EPK-700/1000	Trolley interconnection cable kit
403145	Multi-Video Hybrid Cable 5 m for EPK-700/1000	Trolley interconnection cable kit

Multi-Video Hybrid Cable for EPM 3500

Code	Article	Description
403244	Multi-Video Hybrid Cable 10 m for EPM-3500	Trolley interconnection cable kit
403245	Multi-Video Hybrid Cable 5 m for EPM-3500	Trolley interconnection cable kit

Cable for High Frequency Instruments



Code	Article	Description
83392	AC-ECA	Active Cable for Erbotom: 4 mm \varnothing
83394	AC-VCA	Active Cable for Valleylab: 8 mm \varnothing



Adapter

Olympus-Adapter



Code	Article	Description
83217	AL-OF5	Adapter Olympus OES Scopes to PENTAX Lightsource LH-150(II)/150P(II), LH-150PC, EPM-3000/3300, LX-750P, EPK-1000
83221	AL-OR1	Adapter Olympus rigid scope to PENTAX Lightsource

Code	Article	Description
83245	AL-OL8	Adapter PENTAX FB-X,FCN X to OES CLV-U20, CLK-4



Code	Article	Description
83247	AL-OL9	Adapter PENTAX FB-X, FCN X to OES CLV-U20,CLK-4



Code	Article	Description
83139	AL-OL10	Sleeve, adapter for FNL/FCY/FUR an Olympus CLV-U20

Olympus-Adapter



Code	Article	Description
81358	AL-OL11	Adapter for Olympus Lightsource

Code	Article	Description
83220	AL-OLG6	Adapter PENTAX FS-34P an Olympus Lightsource



Code	Article	Description
83248	AL-OLG9	Adapter PENTAX GI X-series to OES CLV-U20,CLK-4



Code	Article	Description
81359	AL-OLG11	Adapter PENTAX GI V-/W-/P-series to OES CLV-U20,CLK-4

Olympus-Adapter



Code	Article	Description
86639	AL-OLG12	Adapter PENTAX GI X-Series and Ultrasound to OES CLE-145 and CLV-160



Code	Article	Description
86641	AL-OLG13	Adapter PENTAX GI V-/W-/P-series to OES CLE-145



Code	Article	Description
81361	AP-OV1	Adapter for Olympus OES scope to Videomodule PVA-1010; PVK-1000



Olympus-Adapter



Code	Article	Description
83373	AT-OT2	Adapter Lens: Olympus Lecturescope (LS-10) – PENTAX Occular Fiberscopes; 1.0X
83384	AT-OT3	Adapter Lens: Olympus Lecturescope (LS-10) – PENTAX Occular Fiberscopes; 0.8X

Code	Article	Description
86407	AH-OF1	Adapter for Olympus fiberscope connecting PENTAX water bottle OS-H4 to OES





Storz-Adapter



Code	Article	Description
83244	AL-SF1	Adapter Karl Storz Fiberscopes – PENTAX Lightsources LH-150(II), LH-150PC, EPM-3000/3300, LX-750P



Code	Article	Description
83241	AL-SL	Sleeve for PENTAX Endoscopes to Storz Lightsource

Wolf-Adapter



Code	Article	Description
83231	AL-WL	Sleeve for PENTAX Endoscopes to Wolf Lightsource

Adapter for Rigid Endoscopes



Code	Article	Description
83383	AT-ST2	PENTAX Endoscope on Storz Teaching



Code	Article	Description
81363	AP-RV1	Adapter Lens: PENTAX Video Module (PVA-1010/PVK-1000) Rigid Ocular Fiberscopes; 1.0X

Code	Article	Description
83386	AT-UR2	Adapter Lens: PENTAX Observerscope on Rigid Ocular Fiberscopes; 1.0X

Eyepiece Accessory

Code	Article	Description
83043	OF-D1	Eyepiece Cap for PENTAX Endoscopes
83544	OF-D3	Eye shield for PENTAX Fiberscopes; PENTAX Ocular

PENTAX Adapter for Processor, Videomodule, Lightsources

PENTAX Endoscopes



Code	Article	Description
83223	AE-P1	PENTAX Video-Endoscopes on EPM-3000, -3300, -3500

Code	Article	Description
83184	OL-H3	PENTAX Fiber-Endoscopes on LH-150/150P, LH-150II/150PII, LH-150PC, LX-750P



Code	Article	Description
83181	OF-G1	Sleeve for all PENTAX scopes



Code	Article	Description
81362	AP-PV1	Adapter Lens for PENTAX Video Module PVA1010



Accessories for Bedsidescopes

with batterytype



Code	Article	Description
85741	BS-LH1	Lightsource Batterytype for BS-Endoscopes
85742	BS-LH2	Lightsource Batterytype for BS-Endoscopes (brighter type)



Code	Article	Description
85740	BS-LHAC3	Transformer Type Halogen Lightsource 230V for BS-Endoscopes
86080	BS-LHACE	Transformer Type Halogen Lightsource 230V for BS-Endoscopes (brighter type)

Code	Article	Description
85736	BS-LC1	Guide Light Cable for BS-Endoscopes
399230	BS-LB1	Battery 3 V for BS-LH1 (minimum order quantity 5 pcs)
85734	BS-H1	Halogen Lamp Set (2 pcs.) for BS-LH1/BS-LHAC3 bayonet fitting
85735	BS-H2	Halogen Lamp Set (5 pcs.) for BS-LH2/BS-LHACE screw fitting

Code	Article	Description
81953	BS-B1	Cap for BS-Series

Accessories

For halogen light sources



Code	Article	Description
300221	OL-H4	Exchange halogen lamp 150 W for LH-150II/150PII/LH-150PC



Code	Article	Description
300502	OL-H1	Exchange halogen lamp 150W for LH-150/150P/150FP



Autoclavable

Code	Article	Description
84004	OS-H3	Water bottle for SA-P2 (2 pcs.)
86365	OS-H4	Water bottle for PENTAX lightsources and processors
F022-S0002	PLA-BOTTLE OS-H4	Water bottle container for OS-H4

for Processors

Code	Article	Description
86481	OS-A43H	White Balance Adjuster for video processors
84051	SDK-2000	Keyboard for Video Processors EPM-3300/3500
400001	Keyboard	Keyboard incl. cover for EPK-700
400002	Cover	Cover for Keyboard for EPK-700
400003	G844100LCAEU	Keyboard for SAFE-3000 with USB-connector
400004	KEYBOARD	Keyboard for EPK-100p
82753	OS-A39	Foodswitch for PSW-2200
82438	OS-A35	Keyboard for PSW-2200 and EPM-3500
86366	OS-A50	Keyboard for EPK-1000/EPK-100p
87756	OS-A64	Keyboard for SAFE-3000
87161	OS-A70	Keyboard for EPK-i
87083	OS-A61	Food switch for EPK-i (2 pedals)

Luer Lock fitting

Code	Article	Description
83633	OF-B55	T-Adapter for Luer-Lock

Gasinsufflation



Code	Article	Description
83647	OF-G11	CO ₂ -gas adapter for all PENTAX G.I. scopes
82239	OF-B130	CO ₂ Adapter

PENTAX Europe GmbH
LIFE CARE
Julius-Vosseler-Straße 104
22527 Hamburg
Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-Mail: medical@pentax.de

PENTAX U.K. Limited
LIFE CARE
Pentax House
Heron Drive, Langley
Slough SL3 8PN
United Kingdom
Tel.: +44 17 53 / 79 27 23
Fax: +44 17 53 / 79 27 94
E-Mail: lifecare.info@pentax.co.uk

PENTAX France S.A.S.
LIFE CARE
112 quai de Bezons
P. B. 204
95106 Argenteuil
France
Tel.: +33 1 / 30 25 75 75
Fax: +33 1 / 30 25 75 76
E-Mail: medical@pentax.fr

PENTAX Nederland B.V.
LIFE CARE
Lage Mosten 35
4822 NK Breda
Netherlands
Tel.: +31 76 / 5 31 30 31
Fax: +31 76 / 5 31 30 00
E-Mail: medical@pentax.nl

PENTAX Italia S.r.l.
LIFE CARE
Via Dione Cassio, 15
20138 Milano
Italy
Tel.: +39 / 02 50 99 58 1
Fax: +39 / 02 50 99 58 60
E-Mail: marketing.lifecare@pentaxitalia.it

**S.I.M. - Sitemas Integrales
de Medicina, S.A.**
LIFE CARE
Calle Mario Cabré, 13
28030 Madrid
Spain
Tel.: +034-91-3016240
Fax: +034-91-7513115

HOYA Corporation
PENTAX Tokyo office
2-36-9, Maeno-cho
Itabashi-Ku
174-8639 Tokyo
Japan
Tel.: +81 33 / 9 60 51 55
Fax: +81 35 / 3 92 67 24

Your PENTAX sales representative:

Specifications are subject to change without prior notice and without any obligation on the part of the manufacturer.

www.pentaxeuropa.com  PENTAX Corporation · ISO 9001/EN ISO 13485 certified

LCM/01/02/08/300065/02

PENTAX